

# PROVOX<sup>®</sup> ActiValve<sup>®</sup>

Provox<sup>®</sup> Brush

Provox<sup>®</sup> Plug

Provox<sup>®</sup> ActiValve<sup>®</sup> Lubricant

## **Prescriber Information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## **Patents and Trademarks**

Provox® and Provox® ActiValve® are registered trademarks owned by Atos Medical AB, Sweden.

For information about protective rights (e.g. patents), please refer to our web page [www.atosmedical.com/patents](http://www.atosmedical.com/patents).

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore; Fabricante; Fabricante; Tillverkare; Producent; Produsent; Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamintojas; Výrobce; Gyártó; Výrobca; Proizvajalec; Producent; Producător; Proizvodač; Proizvođač; Κατασκευαστής; Производител; Üretici; მწარმოებელი; Ἀρτοποιητής; İstehsalçı; Proizvoditelj; Producers; Pengilang; 製造業者; 제조사; 製造商; 製造商; الجهة المصنعة; היצרן

Date of manufacture; Herstellungsdatum; Productiedatum; Date de fabrication; Data di fabbricazione; Fecha de fabricación; Data de fabrico; Tillverkningsdatum; Fremstillingsdato; Produksjonsdato; Valmistuspäivä; Framleiðsludagsetning; Valmistamise kuupäev; Ražošanas datums; Pagaminimo data; Datum výroby; Gyártás időpontja; Dátum výroby; Datum proizvodnje; Data produkcji; Data fabricației; Datum proizvodnje; Datum proizvodnje; Ημερομηνία κατασκευής; Дата на производство; Üretim tarihi; წარმოების თარიღი; Ἀρτοποιηθέν ἡμερομηνία; İstehsal tarihi; Дата изготовления; Tanggal produksi; Tarikh pembuatan; 製造日; 제조일자; 製造日期; 製造日期; תאריך ייצור; تاريخ التصنيع

Use-by date; Verfallsdatum; Houdbaarheidsdatum; Date limite d'utilisation; Data di scadenza; Fecha de caducidad; Data de validade; Används senast; Användes inden; Brukes innen; Viimeinen käyttöpäivä; Sídasti notkunardagur; Kölblík kuni; Derīguma termiņš; Naudoti iki; Datum expirace; Felhasználhatóság dátuma; Dátum najneskoršej spotreby; Rok uporabnosti; Data ważności; Data expirării; Rok uporabe; Rok upotrebe; Ημερομηνία λήξης; Срок на годност; Son kullanım tarihi; ვარგისობა; Ἀποσυντηρησιμότητας ἡμερομηνία; Son istifadə tarixi; Срок годности; Gunakan sebelum tanggal; Tarikh guna sebelum; 使用期限; 사용 기한; 使用截止日期; 使用截止日期; תאריך אחרון לשימוש; تاريخ انتهاء الصلاحية

Batch code; Chargencode; Partijnummer; Code de lot; Codice lotto; Código de lote; Código do lote; Batch-kod; Batchkode; Batchkode; Eränumero; Lotunúmer; Partii kood; Partijas kods; Partijos kodas; Kód šarže; Sarzs kód; Kód šarže; Številka serije; Kod partii; Codul lotului; Šifra serije; Šifra serije; Кодовый партида; Код на партида; Parti kodu; Код партии; Partiya kodu; Код партии; Kode batch; Kod kelompok; ㄈㄨㄗㄗㄎㄨㄛㄊㄛ; 배치 코드; 批次代碼; 批号; קוד אזור; كود الدفعة

Product reference number; Produktreferenznummer; Referentienummer product; Numéro de référence du produit; Numero riferimento prodotto; Número de referencia del producto; Referência do produto; Referensnummer; Produktreferenznummer; Produktreferenznummer; Artikkelnummer; Tilvísunarnúmer vöru; Toote viitenumber; Produkta atsauces numurs; Gaminio katalogo numeris; Referenční číslo produktu; Termék referenciaszáma; Referenčné číslo produktu; Referenčna številka izdelka; Numer katalogowy produktu; Numărul de referință al produsului; Referentni broj proizvoda; Referentni broj proizvoda; Αριθμός αναφοράς προϊόντος; Референтен номер на продукт; Ürün referans numarası; პროდუქტის საკონტროლო ნომერი; Արտադրանքի տեղեկատու համարը; Məhsula istinad nömrəsi; Справочный номер изделия; Nomor referensi produk; Nombor rujukan produk; 製品参照番号; 제품 참조 번호; 產品參考號; 产品参考编号; מספר סימוכין של המוצר; الرقم المرجعي للمنتج

Do not re-use; Nur zum Einmalgebrauch; Niet hergebruiken; Ne pas réutiliser; Non riutilizzare; No reutilizar; Não reutilizar; Återanvänd ej; Må ikke genanvendes; Skal ikke gjenbrukes; Ei uudelleenkäyttettävä; Má ekki endurnýta; Mitte taaskasutada; Neizmantot atkārtoti; Pakartotinaí nenaudoti; Nepoužívejte opakovaně; Ne használja fel újra; Nepoužívejte opakovane; Ne uporabljajte ponovno; Nie używać ponownie; A nu se reutiliza; Nemojte ponovno upotrebljavati; Ne koristiti ponovo; Μη επαναχρησιμοποιείτε; Не используйте повторно; Tekrar kullanmayın; არ გამოიყენოთ ხელახლა; Մեկանգամյա օգտագործման համար; Yenidən istifadə etməyin; Не подлежит повторному использованию; Jangan gunakan ulang; Jangan guna semula; 再利用不可; 재사용 하지 마십시오; 請勿重複使用; 請勿重复使用; מירעד לשימוש חוזר; يجب عدم إعادة الاستخدام

Non-sterile; Unsteril; Niet steriel; Non stérile; Non sterile; No estéril; Não estéril; Osteril; Ikke-steril; Usteril; Steriloimaton; Ösæft; Mittesteriilne; Nesterils; Nesterilu; Nesterilni; Nem steril; Nesterilné; Nesterilno; Produkt niejałowy; Prodnus nesteril; Nesterilno; Nije sterilno; Μη αποστειρωμένο; Нестерильно; Steril değildir; არასტერილური; Πλ υστηρη; Qeyri-steril; Нестерильно; Tidak steril; Tidak steril; 非滅菌; 무균 아님; 未滅菌; 非无菌; לא מעוקר; غير معقم

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## Provox® ActiValve® Voice Rehabilitation System

It is important that your patient is aware of the consequences of having a tracheostoma. Therefore make sure that the patient is informed and understands the Stomal airway management information below (extract from Provox ActiValve patient manual).

### Stomal airway management

#### Emergency situations

It is important that emergency personnel knows that you are a neck breather. This will help them determining that oxygen and rescue breathing need to be administered to your stoma and not to your mouth or nose.

We recommend that you and your clinician download the guidelines for rescue breathing from

[www.atosmedical.com/rescuebreathing](http://www.atosmedical.com/rescuebreathing)

In the package of your Provox ActiValve you will find a bright yellow Emergency card and Provox ActiValve User Card.

We recommend that you keep both the Emergency Card and Provox ActiValve User Card with you at all times.

#### Elective situations

If you need to undergo a procedure that requires intubation (putting a breathing tube in your windpipe), it is very important that the anesthesiologist and the doctor who is performing the procedure are aware that you are a neck breather and that you are using a voice prosthesis. It is important for them to understand that the voice prosthesis should stay in place. If it is removed, fluids from your mouth, food pipe, or stomach may get into your windpipe.

It is very important that the intubation tube is inserted and removed carefully so that the voice prosthesis does not dislodge or come out.

## 1. Descriptive information

### 1.1 Indications for use

Provox ActiValve is an unsterile indwelling voice prosthesis intended for anterograde insertion in a healed puncture for voice rehabilitation after total laryngectomy. The device is intended for patients who are experiencing early leakage with previous voice prostheses (device life less than 4-8 weeks). The device reduces the need for frequent replacements in a majority of users, but not in all.

### 1.2 CONTRAINDICATIONS

Provox ActiValve is NOT intended:

- for insertion in a freshly made puncture,
- to be in situ during MRI-examination (Magnetic Resonance Imaging), or during Radiation Therapy.

### 1.3 Description of the device

Provox ActiValve Voice Prosthesis is an unsterile single use one-way valve that keeps a TracheoEsophageal (TE)-puncture open for speech, while reducing the risk of fluids and food entering the windpipe. It has a higher opening pressure than standard voice prostheses to prevent inadvertent valve opening when the patient breathes or swallows.

The prosthesis is not a permanent implant and needs periodic replacement. It cannot be sterilized.

The device is available in lengths 4,5, 6, 8, 10, and 12,5 mm and the outer diameter is 7,5 mm (22,5Fr). It comes in three different opening pressures - Light, Strong and XtraStrong (See graph on inner back cover). Magnets in the ring and valve flap determine the opening force. The magnetic force applies only in the closed and nearly closed position of the valve flap. Once the valve flap has opened, the airflow resistance drops. Therefore a higher opening pressure does not mean more air pressure is required during speaking.

The device body is made of medical grade silicone rubber and radio-opaque fluoroplastic. No part of the device is manufactured with natural rubber latex. The nickel coated magnets are completely sealed in epoxy glue with no detectable release of nickel. The device can therefore also be used in patients with hypersensitivity to nickel. Care should however be taken not to physically damage the epoxy seal during removal or other manipulation of the device.

## Provox ActiValve Voice Prosthesis package

The Provox ActiValve comes with the following:

### Clinician items (Fig. 1a)

Provox ActiValve Voice Prosthesis

- (1) Valve flap, fluoroplastic (radio-opaque); silicone hinge.
- (2) Esophageal flange.
- (3) Tracheal flange.
- (4) Valve seat, fluoroplastic (radio-opaque)

- (5) Interacting magnets in valve seat and in valve flap (radio-opaque).
- (6) Safety strap.
- (7) Direction identification knobs.
- (8) Size information.
- (9) Outer diameter: 7.5 mm (22.5 Fr).

Insertion tool (for anterograde insertion)

- (10) Inserter
- (11) Loading Tube

Clinicians Instructions for Use

### **Patient items (shall be given to the patient) (Fig. 1b)**

- 1 bottle (8 ml) Provox ActiValve Lubricant for maintenance,
- 2 reusable Provox Brushes. One for cleaning of the inside of the prosthesis and one for application of Provox ActiValve Lubricant and for inserting the Provox Plug,
- 1 Provox Plug (first aid tool for leakage),
- 2 Provox ActiValve User Cards to be filled in by the clinician,
- 1 Emergency Card,
- 1 Patient's Instructions for Use.

Since the opening forces are not apparent upon simple visual inspection the device comes with **Provox ActiValve User Cards** that also provide important information about the voice prosthesis. Please complete the user cards using a water resistant permanent marking pen, and provide the cards to the patient.

The **Emergency Card** provides important information about emergency airway management. Advise the patient to carry a user card and the emergency card with them at all times in case of emergency, consultation with healthcare providers (e.g. radiology), or during travel. If your patient needs to replace any of the cards, you can order a new one by contacting Atos Medical.

## **1.4 WARNINGS**

Re-use and re-processing may cause cross- contamination and damage to the device, which could cause patient harm.

Dislodgement or extrusion of Provox ActiValve voice prosthesis from the TE puncture and subsequent ingestion, aspiration or tissue damage may occur (See also Section 3).

To reduce the risk of dislodgment/extrusion:

- Select proper prosthesis size.
- If laryngectomy tubes or stoma buttons are used, choose devices that do not exert pressure on the prosthesis or catch onto the tracheal flange during insertion and removal.
- Instruct the patient to consult a physician immediately if there are any signs of tissue edema and /or inflammation /infection.
- Inform the patient to use only original Provox accessories of corresponding size for maintenance and to avoid all other kinds of manipulation of the voice prosthesis.

Do not modify Provox ActiValve voice prosthesis or any of its accessories.

## **1.5 PRECAUTIONS**

- During loading of the prosthesis make sure that the esophageal flange is completely folded forward in the loading tube (Fig. 7). Incorrect loading can cause device damage with immediate valve failure and fluid aspiration.
- Carefully assess any patient on anticoagulant treatment or with bleeding disorders for the risk of bleeding prior to prosthesis placement or replacement.
- Aseptic technique should be employed in the management of the voice prosthesis to avoid introducing contaminants into the stoma, TE-puncture or esophagus.
- For internal lubrication of Provox ActiValve, only use Provox ActiValve Lubricant. Other lubricants may damage the prosthesis irreparably, leading to material destruction and valve malfunction.
- In patients with excessive biofilm formation (micro- organisms and yeasts) the prosthesis, and especially the silicone parts of the prosthesis, may develop heavy biofilm formation. Therefore it is important to instruct the patient to review the Adverse Events section in the manual and to contact their clinician if they experience any of the described events.

# **2. Instructions for use**

**Caution: Do not use the product if the package is damaged or opened.**

## **2.1 Voice prosthesis replacement**

### **Preparation**

**Determine the required length** of the voice prosthesis by using the current (old) prosthesis as measuring device. If there is more than 3 mm / ~0.12 inches space between the tracheal flange and the mucosal wall, a shorter prosthesis should be used.

If the prosthesis fit is too tight, a longer prosthesis should be used. Alternatively, the Provox Measure (see Accessories) can be used to determine the length.

**Determine the required strength.** If the short device life on standard voice prostheses is caused by biofilm, select Provox ActiValve Light. If the shorter device life is caused by biofilm and/or inadvertent valve opening, i.e. the valve flap is opened during swallowing or heavier breathing, select Provox ActiValve Strong.

If an early leakage problem persists, replace the prosthesis with Provox ActiValve that has the next higher opening force. If a higher opening force is perceived to negatively impact the patient's speech, the patient may return to the use of a device with a lower opening pressure, taking shorter device life into account.

### **Loading the voice prosthesis into the insertion tool:**

1. Secure the Safety Strap of the prosthesis by threading it into the slit (Fig. 2) of the Inserter Pin. Fixate the safety strap by pulling it into the narrower part of the slit (Fig. 3). **CAUTION:** Always attach the safety strap of the prosthesis to the Inserter Pin. It could otherwise be lost during insertion causing aspiration or ingestion.
2. Place the prosthesis on top of the Inserter with the tracheal side down (Fig.4).
3. Fold the esophageal flange forward by pinching it together between thumb and index finger. Hold the loading tube with your other hand. Insert the folded esophageal flange of the prosthesis into the slit of the loading tube (Fig. 5).
4. Push the Inserter Pin forward about 1 cm (~0.4 inches) (Fig.6).
5. Press the protruding part of the flange forward with your thumb into the loading tube. The esophageal flange should now be completely folded forward (Figure 7).
6. Push the inserter into the Loading Tube until mark line 1 is level with the edge of the loading tube (Fig.8). Provox ActiValve voice prosthesis is now prepared for placement.

**CAUTION:** If the esophageal flange is not completely folded forward, too much force will be needed to press the prosthesis forward which could cause device damage.

### **Inserting the new prosthesis into the puncture (Anterograde replacement).**

**Caution:** Do not "inject" the prosthesis using the inserter like a syringe.

1. Pull out the current prosthesis with a non-toothed hemostat.
2. Gently insert the Loading Tube in the TE-puncture until the tip is completely inside the puncture.
3. Insert the prosthesis by holding the loading tube stable with one hand and push the inserter forward with your other hand until mark line 2 is level with the rim of the loading tube (Fig. 9).
4. Check the prosthesis placement by holding the loading tube and inserter together and pulling them back through the puncture until the esophageal flange reaches the anterior esophageal wall (Fig. 10).
5. Release the prosthesis by removing the loading tube (the inserter remains in position) (Fig. 11). If the tracheal flange does not unfold automatically, rotate and pull the inserter pin, or use a non-toothed hemostat to pull the tracheal flange and rotate this in position.
6. Check proper position by rotating and lightly pulling the tracheal flange. Thereafter place the voice prosthesis so that the safety strap is pointing downwards.
7. Check valve function by asking the patient to speak (Fig. 12) and by observing the prosthesis while the patient drinks water (Fig. 13).
8. Cut off the safety strap using scissors on the outer edge of the tracheal flange (Fig. 14). Provox ActiValve voice prosthesis is now ready for use.

**Note:** If the insertion is not successful, it can be repeated with the same tools and devices as long as the safety strap has not been cut off. **If the procedure has to be repeated more than twice**, add additional water-soluble lubrication inside the loading tube.

### **Replacement by overshooting**

Overshooting means that during the insertion of the voice prosthesis it is initially pushed entirely into the esophagus. This approach can be useful when inserting a short prosthesis (4.5 mm) or when there is a concern whether the esophageal flange can be entirely delivered into the esophageal lumen using the normal insertion procedure.

**Procedure:** Prepare and insert the valve as described in steps 1-3 above but insert the complete prosthesis into the esophagus, by pushing the Inserter Pin past mark line 2. Release the prosthesis as described in step 5. Leave the safety strap attached to the inserter.

Use a non-toothed hemostat to pull the prosthesis into position. **CAUTION:** Do not pull by the safety strap. Slight rotation might be helpful to fully position the prosthesis. The end of the tracheal flange should point downwards into the trachea and the identification knobs must always face the stoma. Cut the safety strap only after the correct fit is obtained.

### **Retrograde replacement**

**Provox GuideWire** (REF 7215) may be used for retrograde insertion of the voice prosthesis in cases where Provox ActiValve cannot be inserted in an anterograde manner. See Instructions for Use accompanying the Provox GuideWire. Retrograde insertion must be performed by a clinician.

## 2.2 Device lifetime

Biofilm (bacteria and candida) overgrowth will occur in almost all patients. This will eventually cause material degradation, leakage through the prosthesis or other incompetence of the valve. Radiotherapy, salivary content, gastroesophageal reflux, and dietary habits may influence how rapidly biofilm forms. It is therefore not possible to predict the device lifetime.

Laboratory testing shows that in the absence of bacteria and yeasts the device will maintain its structural integrity for 12 months of simulated use. Usage beyond this limit is under the sole discretion of the prescriber.

Indications for replacement include leakage through or (partial) blockage of the valve, degradation of the materials, or medical indications such as problems with the puncture tract (See also Section 3).

## 2.3 Prosthesis maintenance

Instruct the patient to clean the prosthesis at least twice a day and after each meal.

Use the **Provox Brush** and if possible also the **Provox Flush** to clean the prosthesis from food remnants and mucus. This helps to increase the device life.

If excessive candida overgrowth occurs treatment with antifungal medications may be considered as a preventive measure. Laboratory testing shows no negative influence on the function of the voice prosthesis when using the following antifungal medications: Nystatin, Fluconazole and Miconazol. Other drugs directly to or within close proximity of the voice prosthesis should be carefully assessed, since an unknown variety of chemical substances may influence the material.

### Lubrication

Lubricating the inner lumen of the prosthesis helps prevent occasional stickiness of the valve that might otherwise occur, typically after sleep or a prolonged time of not speaking. Therefore instruct the patient that after cleaning the voice prosthesis in the morning the patient should lubricate the inner lumen of Provox ActiValve.

Apply a single drop of Provox ActiValve Lubricant on the second Provox Brush; avoid touching Provox ActiValve Lubricant bottle with the bristles to avoid contamination of the bottle. Brush the inner lumen of Provox ActiValve without touching the stomal mucosa. After use, wipe the brush thoroughly with a gauze to remove any excess oil.

Note that the liquid is thick and that the lubricant bottle releases the lubricant very slowly, one drop at a time. It is important to use only one drop of oil each time and to make sure that any superfluous oil is removed before using the brush (no dripping brush).

Provox ActiValve Lubricant can be used until the expiration date or 8 months after opening the bottle, whichever occurs first. Inform the patient to register the date of opening on the bottle using a water resistant permanent marking pen.

After each use Provox Brush and any other Provox accessory used should be cleaned, and then disinfected at least once a day according to their Instructions for Use.

**Note:** Incidental ingestion or aspiration of small amounts of Provox ActiValve Lubricant is harmless, but may cause coughing. No health problems have been reported regarding repeated aspiration of small amounts of silicone oil, but it is still better to be cautious and avoid it.

## 2.4 Accessories

**CAUTION:** Only use genuine Provox accessories that are intended for use with Provox ActiValve. Other accessories may cause patient harm or product malfunction.

### Accessories to facilitate insertion of the prosthesis (for clinician's use)

- Provox Measure: Indicates the length (prosthesis size) of the TE puncture.
- Provox Dilator: Is used to upsize smaller TE punctures to allow fitting of Provox voice prostheses. The dilator may also be used to temporarily block or stent the TE puncture.
- Provox GuideWire: For placement of a Provox ActiValve, Provox2 or Provox Vega voice prosthesis in a freshly made TE puncture, or for retrograde replacement of a Provox indwelling voice prosthesis.

### Accessories to maintain the device function (for patient's use)

- Provox Brush / Provox Brush XL are used to clean the lumen and valve of the prosthesis.
- Provox Flush allows for flushing the lumen of prosthesis with air or water.
- Provox Plug is a First Aid tool to temporarily block leakage through the prosthesis.

## 2.5 Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

### 3. Adverse events/Trouble shooting information

Indications for replacement of Provox ActiValve voice prosthesis include leakage through the valve, blockage of the prosthesis, bacterial and candidal overgrowth leading to degradation of the materials and /or excessive pressure needed to obtain speech. Other reasons for replacement include medical indications such as problems with the puncture tract.

**Instruct your patient to consult a clinician if any of the following occurs:**

- Any changes in the appearance of the material of the prosthesis or the way it fits in the puncture,
- there is leakage during drinking and cleaning the prosthesis does not help,
- speaking becomes difficult (requires higher effort) and cleaning does not help,
- there are signs of inflammation or tissue change at or near the puncture tract,
- there is bleeding or overgrowth of tissue around the device,
- the patient experiences persistent pain or discomfort in the region of the Provox ActiValve,
- the patient has chronic cough, trouble breathing, or blood in the mucus. These can be signs of a serious health condition, for example rare lipoid pneumonia, or cancer that requires medical attention.

#### Prosthesis related

**Problem:** Leakage through the voice prosthesis. This is a normal development and may be an indication of the end of the device life.

**Solution:** Use a Provox Brush and/or Provox Flush to clean the voice prosthesis. If cleaning of the device does not resolve the leakage, replace the voice prosthesis.

Studies suggest that adding certain food items to the laryngectomized patients' diet may reduce the occurrence of yeasts and bacteria in voice prosthesis biofilms. Food items that have been suggested to have such an effect are *buttermilk*<sub>1,2</sub>, a probiotic milk-drink containing *Lactobacillus Casei Shirota (Yakult®)*<sub>2</sub>, *probiotic bacteria Lactibacillus lactis 53* and *Streptococcus thermophilus*

*B<sub>3,4</sub>* and *Lactobacillus fermentum B54*, *Lactobacillus rhamnosus744*, *Lactobacillus lactis cremoris SK11*, and *Lactobacillus casei Shirota*<sub>4</sub>, and regular caffeinated soft drinks ("Cola")<sub>3</sub>.

**Problem:** Speaking is sometimes not possible (valve does not open).

**Probable cause:** The valve lid sticks to the valve seat, preventing the valve from opening. This becomes more likely the longer the valve has not been in use (e.g. overnight).

**Solution:** Clean the voice prosthesis and lubricate it (see section Cleaning and maintenance). If cleaning and lubricating does not help: Replace the prosthesis.

**Problem:** Leakage through the voice prosthesis immediately after insertion.

**Probable cause:** Immediate leakage after voice prosthesis replacement can occur if the esophageal flange of the prosthesis has not fully deployed into the lumen of the esophagus or if the prosthesis is too short. If the tip of the voice prosthesis inserter has not been placed far enough into the puncture during the insertion, the esophageal flange can unfold in the puncture tract instead of in the esophagus. Immediate leakage after voice prosthesis replacement can also occur if the voice prosthesis has been damaged during the insertion.

**Solution:** Rotate the voice prosthesis to check the fit. Perform an endoscopic examination to determine if the voice prosthesis is in the correct position and of the right length. Correct the position of the device. If the voice prosthesis is too short, replace with a longer device. If the leakage remains after repositioning, pull out the voice prosthesis with the inserter, or a non-toothed hemostat. Inspect the voice prosthesis. If the voice prosthesis has been damaged (e.g. the blue ring has been dislodged or the valve has been torn off), discard the prosthesis, otherwise reload and repeat the insertion procedure as described in section 2.1.

#### Handling related

**Problem:** It is difficult to push the prosthesis forward in the tube.

**Probable cause:** The esophageal flange has not folded completely forward into the loading tube.

**Solution:** Interrupt the insertion and pull the prosthesis back out of the loading tube. Inspect the voice prosthesis. If the voice prosthesis is damaged, discard the prosthesis! If the valve is undamaged reload and repeat the procedure as described in section 2.3.

**Problem:** The prosthesis has been completely inserted into the esophagus (accidental overshooting).

**Solution:** Leave the safety strap attached to the inserter and follow the instructions in section 2.4.

#### Puncture related

**Aspiration of the prosthesis or other components** – Immediate symptoms may include coughing, choking or wheezing.

**Solution:** Partial airway obstruction or complete airway obstruction requires immediate intervention for removal of the object. If the patient can breathe, coughing may remove the foreign body. As with any other foreign body, complications from aspiration may cause obstruction or infection including pneumonia, atelectasis, bronchitis, lung abscess, bronchopulmonary puncture and asthma.

If aspiration of the device is suspected, a CT scan of the lungs should be performed to confirm aspiration and locate the device. If the CT scan confirms aspiration of the device, the device may be retrieved endoscopically using non-toothed grasping forceps. The silicone housing of Provox ActiValve can also be located endoscopically. On a CT scan and during endoscopy, the device may appear as a circular shape with an opening in the middle with an outer diameter of about 14 mm (the flanges of the device), or as a cufflink shape with a shaft length of 4.5, 6, 8, 10 or 12.5 mm, depending on the size of the device. During endoscopy, reflections from the light source on the clear silicone rubber may be seen. Also, in prostheses that have been in situ for some time, white or yellow appearing Candida deposits may be visible on the device.

**Ingestion of the prosthesis** or other components can cause symptoms that largely depend on size, location, degree of obstruction (if any).

**Solution:** The object usually passes spontaneously into the stomach and subsequently through the intestinal tract. If bowel obstruction, bleeding or perforation occurs, or the object fails to pass through the intestinal tract, surgical removal must be considered. If the device has remained in the esophagus it can be removed endoscopically. Spontaneous passage of the device may be awaited for 4-6 days. The patient should be instructed to observe the stools for the ingested device. If the device does not pass spontaneously, or if there are signs of obstruction (fever, vomiting, abdominal pain) a gastroenterologist should be consulted.

The silicone housing of Provox ActiValve can be located and retrieved endoscopically. For details on endoscopic retrieval see section above.

**Hemorrhage/Bleeding** of the puncture – Slight bleeding from the edges of the TE-puncture may occur during replacement of the prosthesis and generally resolves spontaneously. Patients on anti-coagulant therapy, however, should be carefully evaluated for the risk of hemorrhage prior to placement or replacement of the prosthesis.

**Infection and/or edema of the TE-puncture** may increase the length of the puncture tract (e.g during radio therapy).

This may cause the prosthesis to be drawn inward and under the tracheal or esophageal mucosa. Inflammation or overgrowth of the esophageal mucosa may also cause the prosthesis to protrude from the puncture on the tracheal side.

**Solution:** Temporarily replace the prosthesis with a longer one. Treatment with antibiotics with or without corticosteroids may also be considered. If the situation does not improve the prosthesis should be removed. Stenting the puncture with a catheter might be considered.

**Granulation around the TE-puncture.**

**Solution:** Electrical, chemical, or laser cauterization of the area of granulation may be considered.

**Hypertrophic scarring** around the puncture with bulging of the tracheal mucosa over the tracheal flange may occur if the prosthesis is relatively short.

**Solution:** The excess tissue may be removed by using a laser, or a prosthesis with a longer shaft can be inserted.

**Protrusion and subsequent extrusion** of the prosthesis is sometimes observed during infection of the TE-puncture.

**Solution:** Removal of the prosthesis is required to avoid dislodgement into the trachea. Repuncture may be necessary if the puncture has closed.

**Leakage around the prosthesis** – Transient leakage may occur and can improve spontaneously.

**Solution:** The most common reason for persistent leakage around the voice prosthesis is that the device is too long. Replace the prosthesis with a shorter one. If the leakage is indicated by an enlarged puncture, place a Provox XtraFlange over the voice prosthesis' tracheal flange or replace it with a Provox Vega XtraSeal. Temporary removal of the prosthesis, a submucosal purse string suture or insertion of a cuffed tracheal cannula and/or nasogastric feeding tube may also be considered to reduce the size of the puncture. If the leakage persists, surgical closure of the puncture may be necessary.

**Tracheal and/or esophageal tissue damage** may occur if the prosthesis is too short, too long, or is pushed frequently against the esophageal wall by a tracheal cannula, stoma button, or the patient's finger. In patients undergoing radiotherapy this may happen more easily.

**Solution:** If the patient experiences soreness or pain in the area around the voice prosthesis, inspect the tissue around the TE puncture by endoscopy to avoid severe damage.

## 4. Additional information

### 4.1 Date of printing and user assistance information

The version number of this Instruction for Use and manufacturer contact information is located on the back cover.

## 5. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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








**Intended Use**

Provox Brush is a single patient use device intended for cleaning Provox Voice Prosthesis, for insertion of Provox Plug and Provox Vega Plug, for applying lubricant and anti-candida agents and for cleaning of fenestration holes on Provox LaryTubes. The product is intended for use by the patient.

**Before use**

Always inspect the Provox Brush (fig. 1) to make sure that it is not broken or worn out.

Before each use, make sure that

- the protective tip is not cracked or loose.
- the bristles do not appear worn or are loose.
- the wires are not bent or broken.
- the entire handle is not cracked or broken and that the safety wings have not been removed.

If any defects are found, the brush must be discarded.

**Cleaning and disinfection**

Clean the Provox Brush after each use according to fig. 4-7 using drinking water.

Disinfect the Provox Brush (fig. 8) with one of the following methods:

- Ethanol 70% for 10 minutes
- Isopropylalcohol 70% for 10 minutes
- Hydrogenperoxide 3% for 60 minutes

Do not use the device until it is completely dry

**Note:** Always clean and disinfect before use.

**Precaution**

- The device is intended for single patient use only. Reuse in another patient may cause cross contamination.
- Do not clean or disinfect by any other method as this might cause product damage and injury.

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



## 1. Descriptive Information

### 1.1 Intended use

Provox Plug is intended to seal the inner lumen of Provox Voice Prostheses, Provox2 Voice Prostheses and Provox ActiValve and therefore stop leakage of both air and fluids through the voice prosthesis.

### 1.2 Description of the device

Provox Plug is a first-aid tool for temporarily stopping leakage through the voice prosthesis. The device is inserted into the opening of Provox voice prosthesis. The medallion end can be taped to the skin if desired.

### 1.3 Precaution

Provox Plug is intended for single patient use.  
Reuse in another patient may cause cross contamination.

## 2. Instructions for use

### 2.1 Operation instruction

Check device integrity prior to each use by pulling the strap see fig. 1.  
For instructions on how to use the product, please see fig. 2-6.

### 2.2 Cleaning and disinfection

Clean Provox Plug after each use, or if it is changed during use according to Fig. 7-10.  
Fig. 8, Rub the device between your fingers for 15 seconds.

**DO NOT** use any water other than drinking water to clean and rinse the device.  
Disinfect Provox Plug at least once a day if it is removed (Fig. 11) with one of the following methods:

- Ethanol 70% for 10 minutes
  - Isopropylalcohol 70% for 10 minutes
  - Hydrogenperoxide 3% for 60 minutes
- CAUTION: DO NOT** clean or disinfect by any other method as this might cause product damage and injury.

If Provox Plug looks dirty or has air dried in an area with a risk of contamination the device should be both cleaned and disinfected before use.

A risk of contamination could exist if the device has been dropped on the floor or if it has been in contact with a pet, someone with a respiratory infection, or any other gross contamination.

During hospitalization it is important to both clean and disinfect Provox Plug after use but also before use since there is an increased risk of infection and contamination. In a hospital it is better to use sterile water for cleaning and rinsing, rather than drinking water.

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



## **Indications for use**

For use with Provox ActiValve only.

Lubricating the inner lumen of the Provox ActiValve prosthesis helps to prevent sticking of the valve that might otherwise occur e.g. after sleep.

## **Instructions for use**

Once a day, after cleaning in the morning, lubricate the inner lumen of the Provox ActiValve prosthesis with one drop of Provox ActiValve Lubricant using a different Provox Brush. See illustrations.

Do not put more than one drop on the brush. Avoid touching the brush with the ActiValve Lubricant bottle when dispensing a drop. Be patient when squeezing the lubricant bottle. The lubricant bottle releases the thick silicone oil very slowly. If you get Provox ActiValve Lubricant on your skin wash it off with soap and water.

After use clean the brush according to the Instructions for use that accompanies the Provox Brush.

The ActiValve Lubricant can be used until expiration date or 8 months after opening the bottle, whichever occurs first. Register the date of opening on the bottle using a water resistant permanent marking pen.

Note: Avoid touching the stomal mucosa with the brush during application of the Provox ActiValve Lubricant.

Incidental ingestion of small amounts of Provox ActiValve Lubricant is harmless. Incidental aspiration of small amounts of Provox ActiValve Lubricant is also harmless, but may cause coughing.

No health problems have been reported regarding repeated aspiration of small amounts of silicone oil, but it is still better to be cautious and avoid it. Use only one drop of Provox ActiValve Lubricant each time and make sure that any superfluous silicone oil is removed before using the brush (no dripping brush).

## **Disposal**

Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

## **CAUTION:**

- Overuse of the lubricant in quantity or frequency can result in severe coughing, pneumonia or serious breathing problems.
- Do not apply the Provox ActiValve Lubricant directly from the bottle to the prosthesis since this could deposit the liquid in the lungs and cause severe coughing, pneumonia or serious breathing problems.
- Do not enlarge the bottle opening of the Provox ActiValve Lubricant.
- Keep the Provox ActiValve Lubricant out of reach of children and pets that may be harmed by inhaling or swallowing it.

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© Atos Medical AB, Sweden.